

Stainless Steel Diaphragm Valves



**The innovative line of products for
future-oriented technologies in the life science,
biotechnology, cosmetic and food industries**

Take advantage of Georg Fischer's materials and applications know-how for your systems

Quality already begins in the development phase when production concepts and possible risks and consequences are discussed.

That is why we suggest taking advantage of Georg Fischer's materials and applications know-how at an early stage in the design of your system. Complex piping systems are our speciality.





Perfectly adapted valve bodies, diaphragms, bonnets and actuators

The technical data are not binding. They are not warranted characteristics and are subject to change. Please consult our General Conditions of Supply.



The ingenious design of these diaphragm valves guarantees reliable operation even under extremely demanding conditions. Valve bodies, diaphragms, bonnets and actuators are perfectly adapted; the diaphragm or

other single parts are easy to replace. Depending on the application, the bonnet shell is available in stainless steel, cast iron or plastic. Diaphragm valves are ideal for CIP and SIP processes.

Reliable control for demanding processes

Stainless steel diaphragm valve,
Hot Water Sanitization manually operated

Page 6



Stainless steel diaphragm valve, Steam
Sanitization manually operated

Page 7



Stainless steel diaphragm valve, Hot
Water Sanitization Pneumatically actuated

Page 8



Stainless steel diaphragm valve, Steam
Sanitization pneumatically actuated

Page 9-15



Valve bodies in stainless steel

Page 16-17



Multi-port valves in stainless steel

Page 18



Diaphragms

Page 19





Stainless steel diaphragm valve - Hot Water Sanitization

The stainless steel diaphragm valve is a manually operated diaphragm valve with a plastic bonnet. In addition to the compact design and


the position indicator in standard versions, the valve locking option is a valuable extra, useful in all types of diaphragm valve applications.

Product Features


- **Compact, manually operated bonnet in plastic**
- **Light weight**
- **Low torque for easy operation**
- **Sturdy position indicator**
- **Ergonomic hand-wheel**
- **«Clip hand wheel» – easy to remove**

Technical data

Nominal diameter:	DN15–DN80
Valve body:	Forged stainless steel 1.4435/316L or precision casting CF3M/316L
Diaphragms:	Elastomers or PTFE
Bonnet:	Polypropylene, (PPGF)
Hand wheel:	Polypropylene (PPGF/casted stainless steel on request)
Working pressure:	max.10bar
Working temperature:	max.85 °C
Ambient temperature:	max.30°C



Self-draining position



Installation angle
[Depending on dimension]15–45°

Dimensions and Weights

DN mm	D	H	Weight kg	
15	80	90	0.5	
20	80	102	0.9	
25	94	118.5	1.4	
40	117	139	3.0	
50	152	172	4.2	
65		on request		
80		on request		

Stainless steel diaphragm valve Steam Sanitization

Description

The ST 195-MA stainless steel diaphragm valve is a manually operated diaphragm valve with a plastic bonnet. In addition to the compact design and the position indicator in standard versions,

the valve locking option is a valuable extra, useful in all types of diaphragm valve applications.




Technical data

Nominal diameter:	DN15-DN80
Valve body:	Forged stainless steel 1.4435/316L or precision casting CF3M/316L
Diaphragms:	Elastomers or PTFE
Bonnet:	Polypropylene, fiberglass reinforced (PP-GF 30)
Hand wheel:	Polypropylene (PPS/casted stainless steel on request)
Working pressure:	max. 10 bar
Working temperature:	max. 130 °C

Product Features


- **Compact, manually operated bonnet in plastic**
- **Lightweight**
- **Low torque for easy operation**
- **Sturdy position indicator**
- **Ergonomic hand-wheel**
- **«Clip handwheel» - easy to remove**
- **Accessory: locking set to lock the valve**
- **Variable flow direction and installation position**



Dimensions and Weights

DN mm	D	H	Weight kg	
15	80	90	0.5	
20	80	102	0.9	
25	94	118.5	1.4	
40	117	139	3.0	
50	152	172	4.2	
65		on request		
80		on request		

Self-draining position



Installation angle
(Depending on dimension) 15-45°



Stainless steel diaphragm valve - Hot Water Sanitization

This compact diaphragm valve is low-maintenance, has a pneumatic actuator in plastic, and is available in the fail safe to open, fail safe to close and double acting modes of operation. The pneumatic actuators are available in sizes 025 - → 10 bar ← both sides / 028 - → 10 bar at one side. The 028, which is dimensionally

reduced to a minimum, is especially suited to applications where space is limited. A wide range of accessories, such as electrical feedback unit, positioner or stroke limiter, allow optimal adaptation to all types of control tasks.

* ASI interface can be given.

Product Features

- **Compact, pneumatic actuator in plastic**
- **Lightweight**
- **Fail safe to open, fail safe to close or double acting modes of operation**
- **Position indicator is standard**
- **Low maintenance**
- **Wide selection of accessories: electrical feedback unit, stroke limiter and positioner**
- **Variable flow direction and installation position**
- **Autoclavable**

Technical data

Nominal diameter:

DN15–DN80

Valve body:

Forged stainless steel 1.4435/316L or precision casting CF3M/316L

Diaphragms:

Elastomers or PTFE

Actuator housing DN15–DN50:

CoPAGF

Actuator housing DN65/DN80:

PP-GF

Working pressure:

max. 10 bar (025 - → 10 ← / 028 - → 10)

Working temperature:

max. 85°C

Control medium:

Compressed air (oil-free)/neutral, non-aggressive gases, max. 60°C

Max. permissible control pressure:




6 bar for FC (fail safe spring to close) mode,

5 bar for FO (fail safe spring to open) and DA (double acting) modes

Valve design (DN15–DN50)

1. Optical position indicator
2. All-plastic housing are PPGF and housing bottom of stainless steel
3. Preloaded spring sets

4. Lifting spindle assembly of stainless steel
5. CoPA pistons
6. Connection for control air
7. Diaphragm
8. Stainless steel valve body



Stainless steel diaphragm valve Steam Sanitized

This compact diaphragm valve is low- maintenance, has a pneumatic actuator in plastic, and is available in the fail safe to open, fail safe to close and double acting modes of operation. The pneumatic actuators are available in size 025 or 028. The 028, which is dimensionally reduced to a minimum, is especially suited to

applications where space is limited. A wide range of accessories, such as electrical feedback unit, positioner or stroke limiter, allow optimal adaptation to all types of control tasks.

* ASI interface can be given.



Technical data


Nominal diameter:	DN15–DN80
Valve body:	Forged stainless steel 1.4435/316L or precision casting CF3M/316L
Diaphragms:	Elastomers or PTFE
Actuator housing DN15–DN50:	CoPAGF
Actuator housing DN65/DN80:	PP-GF
Working pressure:	max. 10 bar (025 - → 10 ← / 028 - → 10)
Working temperature:	max. 150°C for actuator housing in CoPA max. 130°C for actuator housing in PP
Control medium:	Compressed air (oil-free)/neutral, non-aggressive gases, max. 60°C
Max. permissible control pressure:	6 bar for FC(fail safe spring to close) mode, 5 bar for FO (fail safe spring to open) and DA (double acting) modes

Product Features


- **Compact, pneumatic actuator in plastic**
- **Lightweight**
- **Fail safe to open, fail safe to close or double acting modes of operation**
- **Position indicator is standard**
- **Low maintenance**
- **Wide selection of accessories: electrical feedback unit, stroke limiter and positioner**
- **Variable flow direction and installation position**
- **Autoclavable**

Valve design (DN15–DN50)


1. Optical position indicator
2. All-plastic housing of CoPA or housing top of CoPAGF and housing bottom of stainless steel
3. Preloaded spring sets
4. Lifting spindle assembly of stainless steel
5. CoPA pistons
6. Connection for control air
7. Diaphragm
8. Stainless steel valve body



Fail safe spring to close (FC)



Fail safe spring to open (FO)




Double acting (DA)

Dimensions/Weights, Performance Diagrams, DN15–DN50

Mode: Fail safe spring to close (FC)
 Control pressure: max. 6 bar
 Diaphragms: Elastomers
 Size: 028

Dimensions and Weights



Dimensions for end connections according to ASME BPE connection type

Body with V-band (Hygienic) Clamp Connector (Refer Figure 3)

Connection Size	Flange OD	Pipe ID	Face to Face Length	Max. Valve Height from Pipe Center	Outside Diameter of Actuator
	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)
½" Mini	25.00	9.40	95.5	111.1	68.0
½" Regular	25.00	9.40	109.5	114.6	68.0
¾"	50.40	15.75	109.5	117.7	68.0
1"	50.40	22.10	133.5	161.4	96.0
1½"	50.40	34.80	156.5	234.3	151.0
2"	63.90	47.50	156.5	241.1	151.0
2½"	77.40	60.20	179.5	267.8	180.0
3"	90.90	72.90	179.5	273.7	180.0

Body with V-band (Hygienic) Clamp Connector (Refer Figure 4)

Connection Size	Pipe OD	Pipe ID	Face to Face Length	Max. Valve Height from Pipe Center	Outside Diameter of Actuator
	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)
½" Mini	12.70	9.40	70.0	111.1	68.0
½" Regular	12.70	9.40	84.0	114.6	68.0
¾"	19.05	15.75	84.0	117.7	68.0
1"	25.40	22.10	108.0	161.4	96.0
1½"	38.10	34.80	131.0	234.3	151.0
2"	50.80	47.50	131.0	241.1	151.0
2½"	63.50	60.20	154.0	267.8	180.0
3"	76.20	72.90	154.0	273.7	180.0


Notes:

- For V-band (Hygienic) clamp connections, clamp connector will be welded to the valve body with butt-weld end. Internal finish will be provided as per customer's request.
- Dimensions and tolerances for end connection and face-to-face length are as per ASME BPE.
- All other dimensions are approximate and subject to change without prior notice.

Dimensions/Weights, Performance Diagrams,

Mode: Fail safe spring to open, double acting (FO,DA)
 Control pressure: max. 5 bar
 Diaphragms: all materials
 Size: 025

Dimensions and Weights



Fail safe spring to open/double acting (FO/DA)


DN mm	D mm	D ₁ G	H mm	H ₁ mm	H ₂ mm	H ₃ mm	H ₄ ¹⁾ mm	H ₅ mm	H ₆ ¹⁾ mm	H ₇ mm	H _a mm	L* mm	L ₁ mm	Lift mm	Weight kg
15 (EPDM)	67	1/8"	98	58.5	24	22.5	46	100	155	41.5	73.5		60	8	3.2
15 (PTFE)	94	1/8"	125	68.5	25	22.5	46	100	155	54	86		60	6	3.2
20	94	1/8"	127	70.5	25	22.5	46	100	155	54	86		60	10	3.4
25	94	1/8"	141	84.5	25	22.5	46	100	155	54	86		60	12	3.6
32	116	1/8"	167	97	26	35.5	65	100	160	65	97		60	14	4.8
40	145	1/4"	198	107.5	36	35.5	65	100	170	81	113		60	18	5.6
50	145	1/4"	216	125.5	36	35.5	65	100	170	81	113		60	22	6.6

¹⁾ = with stroke limiter / manual override


*L = for valve bodies see pages 17-18

Control/Working Pressure Diagrams

Fail safe spring to open/double acting (FO/DA) with EPDM diaphragm



with PTFE diaphragm



Maximum working pressure

DN mm	Diaphragms			
	EPDM	PTFE		
	→	→←	→	→←
15	10	10/9	10	10/9
20	10	10/9	10	10/9
25	10	10/9	10	10/9
32	10	10/9	10	10/9
40	10	10/9	10	10/9
50	10	10/9	10	10/9


Working pressure on one side (→)

Working pressure on both sides (→←)

Dimensions/Weights, Performance Diagrams, DN15–DN50

Mode: Fail safe spring to close (FC)
 Control pressure: max. 6 bar
 Diaphragms: all materials
 Size: 025

Dimensions and Weights



Dimensions for end connections according to ASME BPE connection type

Body with V-band (Hygienic) Clamp Connector (Refer Figure 1)

Connection Size	Flange OD	Pipe ID	Face to Face Length	Max. Valve Height from Pipe Center	Outside Diameter of Actuator
	A	B	C	D	E
	(mm)	(mm)	(mm)	(mm)	(mm)
½" Mini	25.00	9.40	63.5	111.1	68.0
½" Regular	25.00	9.40	88.9	114.8	68.0
¾"	25.00	15.75	101.6	117.6	68.0
1"	50.40	22.10	114.3	162.3	96.0
1½"	50.40	34.80	139.7	234.2	151.0
2"	63.90	47.50	158.8	240.4	151.0
2½"	77.40	60.20	193.8	267.7	180.0
3"	90.90	72.90	222.3	273.7	180.0

Body with V-band (Hygienic) Clamp Connector (Refer Figure 2)

Connection Size	Pipe OD	Pipe ID	Face to Face Length	Max. Valve Height from Pipe Center	Outside Diameter of Actuator
	A	B	C	D	E
	(mm)	(mm)	(mm)	(mm)	(mm)
½" Mini	12.70	9.40	63.5	111.1	68.0
½" Regular	12.70	9.40	84.0	114.8	68.0
¾"	19.05	15.75	84.0	117.6	68.0
1"	25.40	22.10	108.0	162.3	96.0
1½"	38.10	34.80	131.0	234.2	151.0
2"	50.80	47.50	131.0	240.4	151.0
2½"	63.50	60.20	154.0	267.7	180.0
3"	76.20	72.90	154.0	273.7	180.0

Notes:


1. Dimensions and tolerances for end connections and face-to-face length are as per ASME BPE.
2. All other dimensions are approximate and subject to change without prior notice.

Dimensions/Weights, Performance Diagrams, DN65-DN80

Mode: Fail safe spring to close, to open, double acting (FC,F0,DA) Control pressure:
max.6 bar for FC mode
max.5 bar for F0 and DA modes

Diaphragms: all materials

Dimensions and Weights



Mode of operation (FC,F0,DA)


DN mm	D mm	D ₁ G	H mm	H ₁ mm	H ₂ mm	H ₃ mm	H ₄ ¹⁾ mm	H ₇ mm	L* mm	Lift mm	Weight kg
65	280	1/4"	344	148	24	46	100	298		25	9.7
80	280	1/4"	348	150	24	46	100	302		35	11.8

¹⁾ = with stroke limiter / manual override


*L = for valve bodies see pages 17-18

Control/Working Pressure Diagrams

Fail safe spring to close (FC)



Fail safe spring to open/double acting (F0/DA)



Maximum working pressure

DN mm	Diaphragms			
	EPDM		PTFE	
	→	→←	→	→←
65	10	10/9	10	7/6
80	10	10/9	10	6/5

Working pressure on one side (→)

Working pressure on both sides (→←)

Accessories for pneumatic actuators, DN15-DN50

Electrical Feedback Unit, ER 52



Type	Type of switch		max. switch.capacity	Code	
ER 52-1	Microswitch	AG,Ni	250V~/6A	199 190 305	
ER 52-2	Microswitch with gold contact	Au	4-30V=/1-100mA	199 190 306	
ER 52-3	Inductive switch	NPN	4.75-30V=/0.1A	199 190 307	
ER 52-4	Inductive switch	PNP	4.75-30V=/0.1A	199 190 308	
ER 52-5	Inductive switch	Namur	8V=	199 190 309	

Complete adaptor kit

The complete adaptor kit must be used if the valve with stroke limiter/manual over ride is combined with the electrical feedback unit ER52.



Dimension	Mode FC,F0,DA Size 028	Mode FC Size 025
DN 15	199 190 387	199 190 387
DN 20	199 190 387	199 190 387
DN 25	199 190 387	199 190 388
DN 32	199 190 388	199 19 0389
DN 40	199 190 389	199 190 389
DN 50	199 190 389	199 190 389

Stroke limiter/Manual override



Dimension	Mode FC,F0,DA Size 028	Mode FC Size 025
DN 15	199 190 381(EPDM) 199 190 382(PTFE)	199 190 382
DN 20	199 190 382	199 190 382
DN 25	199 190 382	199 190 383
DN 32	199 190 383	199 190 384
DN 40	199 190 384	199 190 385
DN 50	199 190 384	199 190 385

Electric pneumatic positioner (for steam sanitisation actuators) Echardt/Sipart



Dimension	Service voltage	Code
EckhardtSRI986	single acting	199 190 348
EckhardtSRI986	double acting	199 190 349
Kit		199 190 350

Electric pneumatic positioner (for hot sanitisation actuators) DSR 500



Discreption	Service voltage	Code
Positioner DSR 500-1	single acting-FC- 24 v/dc	199 190 001
Positioner feedback card	4-20 ma feedback	199 190 575

Accessories for pneumatic actuators, DN65/DN80

Electrical Feedback Unit, ER53

Type	Type of switch		max. switch.capacity	code
ER53-1	Microswitch	AG,Ni	250 V ~ / 10A	199 190 293
ER53-2	Microswitch with gold contact	Au	4-30 V = / 1-100 mA	199 190 297
ER53-3	InductiveswitchwithLED	NPN	9.6-55 V= / 0.2A	199 190 294
ER53-4	Inductive switch with LED	PNP	9.6-55 V= / 0.2A	199 190 295
ER53-5	Inductive switch	Namur	8 V =	199 190 296
ER53-6	Microswitch	EExd	250 V ~ / 5A	199 190 298
ER53-7	Analog signal transmitter		15-30 V = 4-20 mA	199 190 299



Stroke limiter/Manual override

Dimension		Code
d75/DN65	Min./Max.strokelimiter,manualoverride	199190318
d90/DN80		

Pilot valves (for all actuators)

Type PV94/ 95

3/2-way solenoid valve to control single acting pneumatic actuators. Installed either via a base plate with hollow screw directly on the actuator or onto multiple manifolds.

Materials: body polyamide/brass, seals NBR, mode of operation C (fail safe to close).

Electrical connection: cable plug ISO/DIN 43650

Dimension	Service voltage	Code
Version for direct mounting PV -94	230 V, 50-60Hz	199 19 0501
	115 V, 50-60Hz	199 190 500
	24V =	199 190 498
	24 V, 50-60Hz	199 190 499
Version for direct mounting PV -95	230 V, 50 - 60 Hz	199 190 535
	115 V, 50 - 60 Hz	199 190 534
	24 V =	199 190 532
	24 V, 50 - 60 Hz	199 190 533
Multiple manifolds	for 4 valves	199 190 271
	for 6 valves	199 190 272
	for 8 valves	199 190 273
Blank plate		199 190 274
Cable plug with LED and Varistor	230 V, AC/DC	198 806 122
	115 V, AC/DC	198 806 121
	24 V, AC/DC	198 806 123




Type 5470

4/2-way solenoid valve to control double acting pneumatic actuators. Installed with a Namur adaptor plate. For dimensions DN15-50 a Namur adaptor plate must be installed between the actuator and the pilot valve.

Materials: body polyamide/brass, seals NBR, mode of operation G.

Electrical connection: cable plug ISO/DIN 43650

Dimension	Service voltage	Code
Version with Namur plate without throttle screw	230V,AC/DC	199 190 302
	115V,AC/DC	199 190 303
	24V,AC/DC	199 190 304



Adaptor plate NAMUR

Dimension	Mode FC,F0,DA Size 028	Mode FC Size 025
DN 15	199 190 378	199 190 378
DN 20	199 190 378	199 190 378
DN 25	199 190 378	199 19 0378
DN 32	199 190 378	199 190 379
DN 40	199190379	199190379
DN 50	199190379	199190379



Valve bodies

Description

The valve bodies can be assembled to complete valves with all diaphragms, bonnets or pneumatic actuators in our product range. The quality of the surface finish corresponds to

the particular requirements and application. Forged valve bodies can be supplied in all the connection types and dimensions.


Product Features

- **Buttweld ends according to all standards for manual and orbital fusion**
- **Seamless integrated clamp connections, sterile and dairy threads to all standards**
- **Zero static design**
- **Self-draining when positioned accordingly**
- **Ideal for CIP and SIP processes**
- **Practically turbulence free flow**
- **Special designs for individual requirements**

Flow

characteristics

(average values for diaphragm valves)



Technical data

Nominal diameter: DN8–DN80

Version: Stainless steel 316L
cast

forged investment

Δ -ferrite-content: 1.4435 CF8M
 $\leq 0.5\%$ $\leq 1.0\%$

Surface finish:
(average roughness)

Ra = 0.4 μ m Ra = 0.4 μ m

Kv/Cv values (Pipe standard)

Lift in %	kv in l/min., $\Delta p = 1$ bar						
	DN 15	DN 20	DN 25	DN 40	DN 50	DN 65	DN 80
100	70	146	218	684	1156	1571	2533
90	68	140	210	667	1116	1480	2515
80	67	133	201	625	1076	1445	2462
70	64	124	194	610	994	1365	2312
60	59	101	165	545	893	1210	2170
50	47	87	142	457	750	1044	1925
40	35	55	115	345	606	835	1565
30	22	43	65	310	424	625	845
20	17	20	25	180	222	280	401
10	5	7	14	50	64	125	195

Lift in %	Cv in US Gallons/min., $\Delta p = 1$ psi						
	DN 15	DN 20	DN 25	DN 40	DN 50	DN 65	DN 80
100	4.9	10.2	15.3	47.9	81.0	110.0	177.4
90	4.8	9.8	14.7	46.7	78.2	103.6	176.1
80	4.8	9.3	14.1	43.8	75.4	101.2	172.4
70	4.5	8.7	13.6	42.7	69.6	95.6	161.9
60	4.1	7.1	11.6	38.2	62.5	84.7	152.0
50	3.3	6.1	9.9	32.0	52.5	73.1	134.8
40	2.5	3.9	8.1	24.2	42.4	58.5	109.6
30	1.5	3.0	4.6	21.7	29.7	43.8	59.2
20	1.2	1.4	1.8	12.6	15.5	19.6	28.1
10	0.4	0.5	1.0	3.5	4.5	8.8	13.7

Kv/Cv values (Pipe standard)


Lift in %	kv in l/min., $\Delta p = 1$ bar						
	DN 15	DN 20	DN 25	DN 40	DN 50	DN 65	DN 80
100	35	99	180	426	914	1395	2066
90	35	99	173	425	914	1394	2066
80	34	97	155	410	913	1393	2058
70	33	95	145	401	910	1311	2030
60	33	94	140	388	882	1210	1950
50	31	90	105	349	776	1090	1753
40	28	76	95	335	665	872	1445
30	19	65	66	265	525	630	940
20	10	45	48	139	370	250	300
10	3	10	10	22	67	88	142

Lift in %	Cv in US Gallons/min., $\Delta p = 1$ psi						
	DN 15	DN 20	DN 25	DN 40	DN 50	DN 65	DN 80
100	2.5	6.9	12.6	29.8	64.0	97.7	144.7
90	2.5	6.9	12.1	29.8	64.0	97.6	144.7
80	2.4	6.8	10.9	28.7	63.9	97.5	144.1
70	2.3	6.7	10.2	28.1	63.7	91.8	142.2
60	2.3	6.6	9.8	27.2	61.8	84.7	136.6
50	2.2	6.3	7.4	24.4	54.3	76.3	122.8
40	2.0	5.3	6.7	23.5	46.6	61.1	101.2
30	1.3	4.6	4.6	18.6	36.8	44.1	65.8
20	0.7	3.2	3.4	9.7	25.9	17.5	21.0
10	0.2	0.7	0.7	1.5	4.7	6.2	9.9




Dimensions

Buttweld ends



		ASMEBPE
DN	L	øDxs
15	110	12.7x1.65
20	119	19.05x1.65
25	129	25.4x1.65
40	161	38.1x1.65
50	192	50.8x1.65
65	218	63.5x1.65
80	256	76.2x1.65

Clamp ends



			ClampASME BPE fortubeASME BPE	
DN	L	L1	ød	øD
15	108	89	9.40	25.0
20	118	102	15.75	25.0
25	127	114	22.10	50.5
40	159	140	34.80	50.5
50	191	159	47.50	64.0
65	216	194	60.20	77.5
80	254	222	72.90	91.0



Multi-port valves in stainless steel

Description

Multi-port valves are manufactured in various dimensions, angles and positions. Both the available bonnets and pneumatic actuators can be mounted.


Product Features

- Several dimensions, angles and positions possible
- With manually operated bonnet or pneumatic actuator
- Butt weld ends according to all standards for manual and orbital fusion
- Seamless integrated clamp connection, sterile and dairy thread to all standards
- Zero static design
- Self-draining when positioned accordingly
- Ideal for CIP and SIP processes
- Practically turbulence free flow
- Special designs for individual requirements

Diaphragm valve: Zero dead leg valve body

Features:

- Compact design with zero dead leg
- Cast and forged/ bar stock option available
- Forged/ bar stock body made from SS 316L block material – no welded components
- No internal fabrication weld for forged/ bar stock body
- Optimum security against microbial growth
- Option of Sampling/injection points in forged/ bar stock body for high purity systems
- Both mechanical polish and electro polish surface finish options available each with less than 0.40 μm Ra (16AARH)



Diaphragms


Diaphragm qualities

Code	Description	Material	Temperature °C		Shore-A hard kg/cm	Tensile str. kg/cm ² min.
F2	EPDM	Ethylene-Propylene-Diene-Rubber as per FDA	-40	140	68-70	150
70	FPM	Viton	-5	150	74	130
80	Si	Silicon (U.S.P. Class 6)	-40	175	68-70	70
93/30	PTFE/IIR	PTFE/Butyl	-5	150		
93/70	PTFE/FPM	PTFE/Viton	-5	175		
	PTFE/EPDM	PTFE/EPDM	-5	150		
	PTFE/EPDM	PTFE/EPDM«LOW CREEP»	-5	150		
	PTFE/IIR	PTFE/BUTYL«LOW CREEP»	-5	150		



Pressure/Temperature Range (Medium)

Diaphragms



Bonnet/
Actuator

-10°	Manual bonnet, PP-GF (ST 195-MA) Pneumatic actuator, PP-GF, DN65/DN80 (ST 195-DS)	130°
-10°	Pneumatic actuator, CoPa, DN15-DN50 (ST 195-DS)	150°
-10°	Manual bonnet, cast iron (ST 195-G)	175°
-40°	Manual bonnet, stainless steel (ST 195-S)	175°

Head Office

Georg Fischer Piping Systems Pvt.Ltd.

B -9, Kopri Village, Above China Vally Hotel,
Opp. Hiranandani Foundation School,
Powai, Mumbai- 400 076
Maharashtra, India
Contact: 022-4007 2012/2000
Vivek.marathe@Georgfischer.com
Branchoffice@Georgfischer.com

Ahmedabad Office

Georg Fischer Piping Systems Pvt.Ltd.

Regional Office, 906, Sukhsagar Complex,
Ashram Road, Usmanpura,
Ahmedabad - 380 014
Gujarath, India
Contact: 079-27561510

Delhi Office:

Georg Fischer Piping Systems Pvt.Ltd.

732-733, 8 Th Floor, West End Mall,
Janakpuri Dist. Centre,
New Delhi - 110058
Contact: 011-45520494 /95/97

Bengaluru Office:

Georg Fischer Piping Systems Pvt.Ltd.

No-9, 1 St Floor, Opp. Hotel Coral Tree,
R. T. Nagar, Main Road
Bangluru - 32
Contact: 080-41135134

Chennai Office:


Georg Fischer Piping Systems Pvt.

Flat No -10, Door No -12, 3 Rd Floor,
Swathi Complex, Bazullah Road,
T Nagar, Chennai - 600 017
Contact: 044-45502403

Service Centre

Georg Fischer Piping Systems Pvt.Ltd

Gala No. J-1, Jai Matadi Compound,
Ranjan Construction,
Thane-Bhiwandi Road,
Kather Village, Bhiwandi - 421302,
Contact: 091-9819195443



Adding Quality to People's Lives

GFIND/14/CAT/SS/MTL/V0.1

©Georg Fischer Piping Systems Pvt. Ltd

Mumbai, 2014

Printed in India

* These products are for sale in indian territory only



GEORG FISCHER
PIPING SYSTEMS